

Hearing Loss Help Co.

Jodi-Vac Description Sheet

jodi® (Consumer)

The first hearing aid vacuum cleaner designed for the consumer. Actual size; 3"x4"x1 1/2". Weight: just under 2 Lbs. 2lbs. Portable to travel with. 12 Volt DC adds worldwide acceptance.

The filtered syringe assembly is replaceable for maintenance. The syringe has a flow indicator showing when the needle is plugged with debris. A toothbrush is included to clear debris blockages from the needle. A reamer is also included to clear blockages that the brush cannot.

Field tested across North America, has proven effective with all hearing aids.



jodi® -PRO

A portable hearing aid vacuum cleaner for the professional. Actual size; 3"x4"x1 1/2". This pro unit has a 2000 hour clinically rated motor giving the unit a lifetime of use.

Two specially milled needles one for ITE, ITC, and CIC instruments. A second designed for BTE tubing and all other earmolds. The same replaceable filtered syringe provides ease of maintenance and infection control. Syringes and replacement needles are available in economical packages. The syringe has a built in flow indicator to make it easy to tell when the needle is plugged.

A brush is molded into the box for dislodging debris from the needle, and a reamer is included for the tough-to-dislodge blockages. Low voltage 12 volt DC adds to the portability.



jodi® -XL (New Improved)

29 PSI air pressure—blower feature. Vacuum sucks 21 inches of mercury.

A clinical vacuum desiccator for the office, clinic, or hospital. Two vacuum pumps power this unique unit.

The desiccator has a heated plate and a desiccant for maximum moisture absorption. The motor turns off when it reaches maximum vacuum pressure until needed to keep the desiccator at vacuum. The second pump powers the needle vac which is identical to the pro unit.

Small in size, powerful and durable enough to last a lifetime, this is the ultimate hearing aid cleaning tool. 15 volt DC with a dual 110-230 transformer. This unit can be used anywhere in the world.

Dimensions: 9.25 "L x 6.25"W x 5"H (with lid)



To order, go to: www.Hearing-Loss-Help-Co.com.

Or Email: HearingLossHelpCo@Hearing-Loss-Help-Co.com.

Write: Hearing Loss Help Co, 75 Townsend St., Worcester, MA 01609-1979.

PH: (508) 755-1120. FAX: 508 546 8000